


Reagent consumption chart



FOCUSTRONIC

SETTING NEW STANDARDS

Comment:

- A measure cycle consumes 6-9 ml of reagent.
- If you have a wash out value (black dots) the reagent consumption behind that value is twice the normal for that measurement (12-18 ml).
- The chart below assumes that you have a minority of such was out measurements (the normal case), and an average dKH around 8dKH. If lower dKH reagent lasts some longer.

Measure interval	4 l conc bottle last (after diluted to ready to use)
2 h	210 days
4 h	420 days
6 h	630 days
8 h	840 days